Luxatemp Fluorescence

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Product description

Luxatemp Fluorescence is a self-curing composite for the fabrication of temporary crowns, bridges, inlays, onlays, partial crowns and veneers. The 2-component material is based on multi-functional methacrylates and it is automatically mixed. The additional fluorescence makes it possible to create particularly natural temporaries of high esthetic quality.

Intended purpose

Fabrication of temporaries and long-term temporaries: crowns, bridges, inlays, onlays, partial crowns and veneers

Limitations of use

Do not use the material for permanent restorations.

Indication

Restoration of oral function and esthetics:

- Major tooth defects that cannot be restored by direct restoration, e.g. extensive carious lesions, traumatic injuries or extensive dental wear
- Endodontically treated teeth that are susceptible to fractures
- Tooth loss
- Dental misalignments and deformities, e.g. individual or multiple diastemas and internal discoloration
- Cases for which it is not possible to predict the functional or esthetic outcomes of direct restorations
- Occlusal modification

Contraindication

Do not use the material where there are known allergies to any of the components or contact allergies.

Patient target group

Persons receiving treatment as part of a dental procedure.

Intended user

Dentist, dental assistant staff and dental technicians

Notes for use

- With the Smartmix version, the material is applied directly without a dispensing device.
- With the Automix version, the material is applied using a 10:1 dispenser.
- When using an Automix Cartridge or a Smartmix Syringe for the first time, squeeze out a pea-sized amount of material and discard it.

Recommended use

- Before preparing the crown or bridge or prior to extraction of a tooth, take an impression using alginate, silicone or a thermoplastic impression material.

Alternatively, the temporary can be made using a laboratory-fabricated vacuum-formed stent or, in the case of single crowns, with a strip crown. Adhere precisely to the relevant manufacturer's instructions!

Note: In a silicone impression, undercuts should be compensated for and, if necessary, relief channels should be cut. To create a connection with the abutment teeth in the case of a gap in the molar region, cut a groove into the casting.

- Dry the prepared teeth and compensate for any undercuts present in the preparation. Then lightly lubricate the cores and surrounding tissue plus any composite core build-ups using petroleum jelly or a similar separating medium.
- Select the desired Luxatemp Fluorescence shade and insert the relevant cartridge into the Automix Dispenser or use the relevant Smartmix Syringe (see "Using the Automix Cartridge"/see "Using the Smartmix Syringe").

Note: To avoid bubbles, it is important that the end of the mixing tip always remains immersed in the material during application.

- Dispense Luxatemp Fluorescence into the impression or vacuum-formed stent. Then dispense the material onto the occlusal surfaces of the impression or vacuum-formed stent and fill up to the gingival areas, only overbuilding slightly.
- Place the impression or vacuum-formed stent filled with Luxatemp Fluorescence on the prepared teeth 45 seconds at the latest after mixing has begun and hold in place with moderate pressure.

Note: The setting process must be monitored intraorally (e.g. using a probe), as the temperature in the mouth has a critical effect on curing, and the provisional can only be removed during its elastic phase.

- The provisional must be removed from the mouth during the elastic phase, i.e. during the period from 2 to 3 min after mixing has begun.

Note: To further improve the accuracy of large temporaries, they can be placed back into the impression for final curing after removal from the mouth.

- Before finishing, remove the oxidation prevention layer from the surface of the provisional (e.g. with alcohol).
- From about 6 to 7 min after mixing has begun, the provisional can be finished using cross-cut grinders, flexible discs, etc.
- Check occlusion and adjust if necessary using suitable instruments.
- If necessary, coat the surface of the provisional with a light-curing high-gloss glaze (e.g. Luxatemp Glaze & Bond) or manually polish it with a rubber polisher or small goat's hair brush.
- Clean the core, dry it and thoroughly remove the separating medium.
- Attach the provisional.

Note: Any temporary cement can be used for fixing the provisional. DMG recommends the use of automatically mixing luting cement. Cements that contain eugenol may prevent the polymerization of resin-based luting cements during the definitive cementation.

- Check the occlusion and adjust accordingly.

Customizing and repairing a provisional

Note: When customizing and repairing a freshly-made provisional, it does not need to be roughened and a bonding agent is not required (continue with step 2).

- Repairing a provisional that is being worn: Roughen the provisional slightly using a grinder or a sandblaster at the appropriate points. Afterwards apply a bonding agent (e.g. Luxatemp Glaze & Bond) in accordance with the manufacturer's instructions.
- Apply Luxatemp Fluorescence to all surfaces.
- Fix parts for 3 min.
- After 6 or 7 min, the appropriate points can be finished as usual.

Note: To speed up the polymerization, the processed provisional can be placed in warm water at 50 °C (120 °F) for a few minutes. Alternatively, the provisional can also be customized or repaired using the light-curing, easy-flowing composite LuxaFlow. Follow the corresponding instructions for use!

Timing

Place into the mouth	0:00–0:45 min
Remove from the mouth	2:00–3:00 min
End of curing	6:00–7:00 min
Processing the provisional	

Note: The indicated times apply at a room temperature of 23 °C/73 °F and a normal relative humidity of 50 %. Higher temperatures shorten, lower temperatures lengthen these times.

Using the Automix Cartridge (for image, see "Handling")

Use of the Automix Cartridge

- Push up the lever [A] at the back of the Automix Dispenser and pull back the slide [B] all the way.
- Lift up the plastic lever [C], insert the Automix Cartridge and lock it using the plastic lever [C].

Note: Make sure that the notches on the Automix Cartridge line up with those on the Automix Dispenser.

Attaching a mixing tip

- Remove the cartridge cap [D] or used mixing tip after turning it 90 ° counterclockwise, and discard it.
- Attach a new mixing tip [E].

Note: Make sure that the notches on the mixing tip line up with those on the cartridge.

- Twist the mixing tip clockwise by 90 ° to lock it.

Applying the material

- When using an Automix Cartridge for the first time, squeeze out a pea-sized amount of material and discard it.
- Mix the material in the mixing tip using the lever [F] on the Automix Dispenser and apply directly.

Note: After the application, leave the used mixing tip on the Automix Cartridge as a seal!

Removal of the Automix Cartridge

1. Push up the lever [A] at the back of the Automix Dispenser and pull back the slide [B] all the way.
2. Lift up the plastic lever [C] and remove the Automix Cartridge.

Use of the Smartmix Syringe (for image, see “Handling”)

1. Before attaching the mixing tip, remove the cap [A] or used mixing tip after turning it 90 ° counterclockwise, and discard it.
2. Attach a new mixing tip [B].

Note: Make sure that the notches on the mixing tip line up with those on the Smartmix Syringe [C].

3. Twist the mixing tip clockwise by 90 ° to lock it.
4. When using a Smartmix Syringe for the first time, squeeze out a pea-sized amount of material and discard it.
5. The material is mixed in the mixing tip when it is expressed and it can be applied directly.

Note: After application, leave the used mixing tip on the Smartmix Syringe as a seal!

Disinfection of the Automix Dispenser

- The Automix Dispenser can be reused after disinfection, although it should be replaced if there are visible signs of wear.
- Disinfection can be carried out using a commercially available dipping disinfectant. DMG recommends using only disinfectant solutions listed by the RKI (Robert Koch Institute). The Automix Dispenser is also autoclavable.

Interactions

- Residues of methacrylate-based materials may alter the setting behavior of the silicone impression and bite registration materials.
- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area.

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

Warnings/precautions

- Keep out of the reach of children!
- For dental use only!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Reuse: Standard surface disinfection as required. No specific maintenance required.
- Follow the manufacturer’s instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.
- When finishing the provisional, sanding dust can get into the eyes and respiratory tract.
- Wear a face mask and protective goggles.
- Protect against grinding dust by proper ventilation – do not inhale the dust.

Product performance characteristics

Compressive strength ≥ 250 MPa, flexural strength ≥ 75 MPa.

Storage/disposal

- Store in a dry place at 15 to 25 °C/59 to 77 °F! Do not store in the refrigerator!
- Protect against exposure to direct sunlight!
- After the application, leave the used mixing tip on the Automix Cartridge as a seal!
- Do not use after the expiration date!
- Disposal must comply with national regulations.

Composition

Dental glass, EBPADMA, aliph. polyester diurethanes, unsaturated polyester resin, SiO₂, additives.

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