Splint PMMA BioStar

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Indications

Splint PMMA BioStar are transparent dental milling blanks out of PMMA for the manufacture of splints, therapeutic splints, bite regulators and drilling templates for long-term application in the oral cavity up to 12 months.

Contraindications

- intolerance to the contained components
- bridges with two molar pontics
- too little connector cross sections
- model casting with clasps, bars, narrow sublingual bars and transversal connectors

Safety information

Please pay attention to the information in the instruction manual in its current version. Avoid the inhalation of dust particles during processing. Wear protection gloves, safety glasses and a dust mask to prevent irritations of eyes, skin and respiratory system.

Handling and storage

Only use clean, defect-free material out of undamaged packaging and make sure that the blanks are always stored in their original container in a dry and cool place, protected from light. Avoid contact with fluids.

Milling and further processing:

Tungsten carbide burs like the one you use for soft reline materials are suitable for rough adaption. For the finalization common cross-toothed millers can be used. Polishing should be carried out with regard to a low contact pressure, otherwise the resulting heat might lead to smearing or melting.

Only use one-edged milling tools in the milling machine or wet grinding procedures!

Cleaning

- ultrasonic bath (max. 40 °C) with mild cleaning agent
- Do not use the following methods for cleaning:
- steam blasting
 - autoclave method (e.g. sterilization / disinfection)
 - organic solvents and acids (e.g. acetone, hydrochloric acid)
 - highly concentrated alcohols
 - chemical denture cleaners
 - abrasive cleaning agents

Disposal

Remains of PMMA blanks can be disposed of in the regular household garbage.

We reserve the right to make changes because of the continuous development of our products. Please find the current version of the instruction manual on our website: www.hinrichs-dental.de This version replaces all previous versions.

Explanation of symbols

Manufacturer

Consult instructions for use

Batch code

Store dry

Caution: U.S. Federal law restricts this device to sale by or on the order of a dentist.

Content (Quantity)

Date of manufacture

Use by

Catalogue No.

Protect from direct sun light

Date of information: 11/2020

Printing date: 19.11.2020

Hersteller / Manufacturer:
Dental Direkt GmbH
Industriezentrum 106 – 108
32139 Spempe, Germany

REV: 11/2020

Distributor:
ERNST HIRCHS Dental GmbH
DE-38644 Goslar · Borsigstr. 1

Tel.: +49 (0) 53 21 / 50 624
Fax: +49 (0) 53 21 / 50 881

info@hinrichs-dental.de
www.hinrichs-dental.de