
IPS Ivocolor®

Combination table for staining the IPS veneering materials and the all-ceramic materials of the IPS e.max® system*

Kombinationstabelle für die Bemalung von IPS Verblendkeramiken und Vollkeramiken des IPS e.max® Systems*



* Other shade combinations apply when using IPS e.max ZirCAD LT discs in the LT group shades.
Bei der Verwendung von IPS e.max ZirCAD LT Discs in den LT-Gruppenfarben gelten andere Farbkombinationen.

see instructions

ivoclar vivadent® technical

IPS Ivocolor®

Combination table for IPS Ivocolor® components

Kombinationstabelle für IPS Ivocolor® Komponenten

	Initial mixing Erstmaliges Anmischen			Refreshing of dried-out materials Auffrischung von eingetrockneten Massen		
Components Komponenten	**IPS Ivocolor Mixing Liquid allround** 	**IPS Ivocolor Mixing Liquid longlife** 	**IPS Ivocolor Essence Fluid** 	**IPS Ivocolor Mixing Liquid allround** 	**IPS Ivocolor Mixing Liquid longlife** 	**IPS Ivocolor Essence Fluid** 
 **IPS Ivocolor Essence**	✓	✓	✓	✓	✓	✗
 **IPS Ivocolor Shade**	✓	✓	✗	✓	✓	✗
 **IPS Ivocolor Glaze Powder/ FLUO**	✓	✓	✗	✓	✓	✗
 **IPS Ivocolor Glaze Paste/ FLUO**	✓	✓	✗	✓	✓	✗

CE 0123

Rx only

Manufacturer:
Ivoclar Vivadent AG, 9494 Schaan/Liechtenstein
www.ivoclarvivadent.com

see instructions

Date information prepared: 2019-10, Rev. 1
668431/d-e

ivoclar vivadent: technical